

## LAP3 Rabbit mAb

Catalog No: #57145

Package Size: #57145-1 50ul #57145-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	LAP3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	LAP3 Antibody detects endogenous levels of total LAP3
Immunogen Description	A synthesized peptide derived from human LAP3
Other Names	LAP; Lap3; LAPEP; Leucine aminopeptidase 3; PEPS; Peptidase S;
Accession No.	P28838
Calculated MW	52kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200

## Images



Western blot analysis of LAP3 expression in 293T cell lysate.

## Background

Presumably involved in the processing and regular turnover of intracellular proteins. Catalyzes the removal of unsubstituted N-terminal amino acids from various peptides.

---

Note: This product is for in vitro research use only