

ROR1 Antibody FITC Conjugated

Catalog No: #C00580F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ROR1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 300-350 937 derived from human ROR1
Conjugates	FITC
Target Name	ROR1
Other Names	NTRKR1; dJ537F10.1; Tyrosine-protein kinase transmembrane receptor ROR1; Neurotrophic tyrosine kinase, receptor-related 1; ROR1
Accession No.	Swiss-Prot#Q01973NCBI Gene ID4919
Uniprot	Q01973
GeneID	4919;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

Tyrosine-protein kinase receptor whose role is not yet clear. $\text{ATP} + \text{a [protein]-L-tyrosine} = \text{ADP} + \text{a [protein]-L-tyrosine phosphate}$.

Note: This product is for in vitro research use only