

DRD2 Antibody FITC Conjugated

Catalog No: #C00586F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DRD2 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 214-260 443 derived from human DRD2
Conjugates	FITC
Target Name	DRD2
Other Names	D2R; D2DR; D(2) dopamine receptor; Dopamine D2 receptor; DRD2
Accession No.	Swiss-Prot#P14416NCBI Gene ID1813
Uniprot	P14416
GeneID	1813;
Excitation Emission	494nm 518nm
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Dopamine receptor whose activity is mediated by G proteins which inhibit adenylyl cyclase.

Note: This product is for in vitro research use only