

Desmoglein 2 Antibody FITC Conjugated

Catalog No: #C00623F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Desmoglein 2 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 60-110 1118 derived from human Desmoglein 2
Conjugates	FITC
Target Name	Desmoglein 2
Other Names	HDGC; CDHF5; ARVC10; ARVD10; CMD1BB; Desmoglein-2; Cadherin family member 5; DSG2
Accession No.	Swiss-Prot#Q14126NCBI Gene ID1829
Uniprot	Q14126
GeneID	1829;
Excitation Emission	494nm 518nm
Cell Localization	Cytoplasm, Nucleus
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

Component of intercellular desmosome junctions. Involved in the interaction of plaque proteins and intermediate filaments mediating cell-cell adhesion.

Note: This product is for in vitro research use only