

Recombinant Human Sushi repeat-containing protein SRPX2(SRPX2), Escherichia Coli

Catalog No: #AP89654

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Sushi repeat-containing protein SRPX2(SRPX2), Escherichia Coli
Brief Description	Recombinant protein
Host Species	E.coli
Purification	>85% (SDS-PAGE)
Immunogen Description	Protein Length: Full Length of Mature Protein Expression Region: 24-465
Target Name	His-tagged
Other Names	BPP; CBPS; PMGX; RESDX; SRPUL; SRPX 2; SRPX2;
Accession No.	Swiss Prot:O60687
Uniprot	O60687
Target Species	Homo sapiens (Human)
Tag Info	The tag type will be determined during production process
Target Sequence	Protein Length: Full Length of Mature Protein Expression Region: 24-465 TWYAGSG YYPDESYNEV YAEVVPQAPA LDYRVPRWCY TLNIQDGEAT CYSPKGGNYH SSLGTRCELS CDRGFRLIGR RSVQCLPSRR WSGTAYCRQM RCHALPFITS GTYTCTNGVL LDSRCDYSCS SGYHLEGDRS RICMEDGRWS GGEPVCVDID PPKIRCPHSR EKMAEPEKLT ARVYWDPPPLV KDSADGTITR VTLRGPEPGS HFPEGEHVIR YTAYDRAYNR ASCKFIVKVQ VRRCP TLKPP QHGYLTCTSA GDNYGATCEY HCDGGYDRQG TPSRVCQSSR QWSGSPICA PMKINNVNNS AAGLLDQFYE KQRLLIISAP DPSNRYKMQ ISMLQQSTCG LDLRHVTIIE LVGQPPQEVG RIREQQLSAN IIEELRQFQR LTRSYFNMVL IDKQGIDRDR YMEPVTPEEI FTFIDYLLS NQELTQRREQ RDICE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°/-80°. The shelf life of lyophilized form is 12 months at -20°/-80°

Note: This product is for in vitro research use only