

PTGR2 Polyclonal Antibody

Catalog No: #47347

Package Size: #47347-1 50ul #47347-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PTGR2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Rat;Mouse
Specificity	PTGR2 Polyclonal Antibody detects endogenous levels of PTGR2 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human ZADH1.
Target Name	PTGR2
Uniprot	Q8N8N7
GeneID	145482
Calculated MW	35kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. ELISA: 1/10000.. IF 1:50-200

Images



Western Blot analysis of K562 cells using PTGR2 Polyclonal Antibody



Immunohistochemical analysis of paraffin-embedded human tonsil.

Note: This product is for in vitro research use only