

Human Ubiquitin thioesterase OTU1 ELISA Kit

Catalog No: #EK21121

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human Ubiquitin thioesterase OTU1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.14ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only