

SLC25A20 rabbit pAb

Catalog No: #48787



Package Size: #48787-1 50ul #48787-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SLC25A20 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of SLC25A20 at Human, Mouse,Rat
Immunogen Description	Synthesized peptide derived from human SLC25A20
Target Name	SLC25A20
Uniprot	O43772
GeneID	788
Calculated MW	33kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.507% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000

Note: This product is for in vitro research use only