

CD134/OX40L Receptor Rabbit mAb

Catalog No: #57237

Package Size: #57237-1 50ul #57237-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CD134/OX40L Receptor Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP
Species Reactivity	Mouse
Immunogen Description	Recombinant protein of mouse CD134
Other Names	OX40 antigen; OX40L receptor
Accession No.	Swiss-Prot:P47741 Gene ID:22163
Uniprot	P47741
GeneID	22163;
Calculated MW	30 kDa
SDS-PAGE MW	45 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000IHC: 1/50-1/100

Images



Western blot analysis of CD134 in mouse spleen lysates using CD134/OX40L Receptor antibody.

Note: This product is for in vitro research use only