

HSD17B1 Rabbit mAb

Catalog No: #57276



Package Size: #57276-1 50ul #57276-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	HSD17B1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP ICC
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human HSD17B1
Other Names	hydroxysteroid (17-beta) dehydrogenase 1; HSD17; EDHB17; EDH17B2; SDR28C1
Accession No.	Swiss-Prot:P14061 Gene ID:3292
Uniprot	P14061
GeneID	3292;
Calculated MW	35 kDa
SDS-PAGE MW	35 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP: 1/20WB: 1/500-1/1000IHC: 1/50-1/100ICC: 1/50-1/200

Images



Western blot analysis of HSD17B1 in Human placenta lysates using HSD17B1 antibody.

Note: This product is for in vitro research use only