

N4BP1 Rabbit mAb

Catalog No: #57295



Package Size: #57295-1 50ul #57295-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	N4BP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human
Immunogen Description	Recombinant protein of mouse N4BP1
Other Names	C81621; AI481586
Accession No.	Swiss-Prot:Q6A037 Gene ID:80750
Uniprot	Q6A037
GeneID	80750;
Calculated MW	99 kDa
SDS-PAGE MW	99 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of N4BP1 in K562, HeLa lysates using N4BP1 antibody.

Note: This product is for in vitro research use only