

Prostaglandin D Synthase Rabbit mAb

Catalog No: #57323



Package Size: #57323-1 50ul #57323-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Prostaglandin D Synthase Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC ICC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human Prostaglandin D Synthase
Other Names	PDS; PGD2; PGDS; LPGDS; PGDS2; L-PGDS
Accession No.	Swiss-Prot:P41222 Gene ID:5730
Uniprot	P41222
GeneID	5730;
Calculated MW	21 kDa
SDS-PAGE MW	26 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB: 1/500-1/1000IHC: 1/50-1/100ICC: 1/50-1/200



Images



Western blot analysis of Prostaglandin D Synthase in K562, C6, 3T3, HeLa lysates using Prostaglandin D Synthase antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using Prostaglandin D Synthase antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of Prostaglandin D Synthase (green) in K562 using Prostaglandin D Synthase antibody, and DAPI (blue).

Note: This product is for in vitro research use only