

LLGL1 Rabbit mAb

Catalog No: #57381



Package Size: #57381-1 50ul #57381-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	LLGL1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human
Immunogen Description	Recombinant protein of human LLGL1
Other Names	DLG4; HUGL; LLGL; Lgl1; Mgl1; HUGL1; HUGL-1
Accession No.	Swiss-Prot:Q15334 Gene ID:3996
Uniprot	Q15334
GeneID	3996;
Calculated MW	115 kDa
SDS-PAGE MW	115 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of LLGL1 in A549, U2OS lysates using LLGL1 antibody.

Note: This product is for in vitro research use only