

Sclerostin Antibody HRP Conjugated

Catalog No: #C00644H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sclerostin Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F ICC
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 150-200 213 derived from human Sclerostin
Conjugates	HRP
Target Name	Sclerostin
Other Names	CDD; VBCH; SOST1; Sclerostin; SOST; UNQ2976 PRO7455 PRO7476
Accession No.	Swiss-Prot#Q9BQB4NCBI Gene ID50964
Uniprot	Q9BQB4
GeneID	50964;
Excitation Emission	N A
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Background

Negative regulator of bone growth that acts through inhibition of Wnt signaling and bone formation.

Note: This product is for in vitro research use only