

Recombinant Human B-cell Lymphoma-extra Large

Catalog No: #AP60468




Package Size: #AP60468-1 10ug #AP60468-2 100ug #AP60468-3 500ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human B-cell Lymphoma-extra Large
Host Species	Escherichia coli.
Purification	> 97 % by SDS-PAGE and HPLC analyses.
Uniprot	Q07817
GeneID	598
Calculated MW	Approximately 23.7 kDa, a single non-glycosylated polypeptide chain containing 211 amino acids.
Target Sequence	SQSNRELVVD FLSYKLSQKG YSWSQFSDVE ENRTEAPEGT ESEMETPSAI NGNPSWHLAD SPAVNGATGH SSSLDAREVI PMAAVKQALR EAGDEFELRY RRAFSDLTSQ LHITPGTAYQ SFEQVVNELF RDGVNWGRIV AFFSFGGALC VESVDKEMQV LVSRIAAWMA TYLNDHLEPW IQENGGWDTF VELYGNNAAA ESRKGQERFN R
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in 20 mM Tris-HCl, pH 8.0, 5 % Trehalose.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Images



Note: This product is for in vitro research use only