

# CD5 Antibody FITC Conjugated

Catalog No: #C00654F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |   |
|-----------------------|---|
| Product Name          | CD5 Antibody FITC Conjugated  |
| Host Species          | Rabbit  |
| Clonality             | Polyclonal  |
| Isotype               | IgG   |
| Purification          | Purified by Protein A.  |
| Applications          | ICC IF  |
| Species Reactivity    | Hu Ms Rt  |
| Immunogen Description | KLH conjugated synthetic peptide aa 250-300 495 derived from human CD5          |
| Conjugates            | FITC  |
| Target Name           | CD5   |
| Other Names           | T1; LEU1; T-cell surface glycoprotein CD5; Lymphocyte antigen T1 Leu-1; CD5     |
| Accession No.         | Swiss-Prot#P06127NCBI Gene ID921  |
| Uniprot               | P06127  |
| GeneID                | 921;  |
| Excitation Emission   | 494nm 518nm   |
| Cell Localization     | Extracellular   |
| Concentration         | 1mg ml  |
| Formulation           | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.                |
| Storage               | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

## Application Details

ICC=1:50-200 IF=1:50-200

## Background

May act as a receptor in regulating T-cell proliferation.

Note: This product is for in vitro research use only