

Flotillin 1 Rabbit mAb

Catalog No: #58682

Package Size: #58682-1 50ul #58682-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Flotillin 1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human Mouse Rat
Specificity	Flotillin 1 Antibody detects endogenous levels of total Flotillin 1
Immunogen Description	A synthesized peptide derived from human Flotillin 1
Other Names	Flotillin1 ; Flotillin 1; Reggie 2; FLOT1; Integral membrane component of caveolae;
Accession No.	Uniprot:O75955
Uniprot	O75955
Calculated MW	48kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50


Images



Western blot analysis of Flotillin 1 expression in (1) HeLa cell lysate; (2) K562 cell lysate.



Immunohistochemical analysis of paraffin-embedded human uterus, using Flotillin 1 Antibody .



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

Product Description

May act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles.

Background

May act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles.

Note: This product is for in vitro research use only