

Cyclophilin B Rabbit mAb

Catalog No: #58767

Package Size: #58767-1 50ul #58767-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Cyclophilin B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Cyclophilin B Antibody detects endogenous levels of total Cyclophilin B
Immunogen Description	A synthesized peptide derived from human Cyclophilin B
Other Names	Cphn2; Cyclophilin B; OI9; CYPB; PPIase B; PPIB; SCYLP; Rotamase B;
Accession No.	Uniprot:P23284
Uniprot	P23284
Calculated MW	18kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB 1:500~1:2000

Images



Western blot analysis of Cyclophilin B expression in HepG2 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

Product Description

PPlases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides.

Background

PPlases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides.

Note: This product is for in vitro research use only