

IRF6 Rabbit mAb

Catalog No: #59072



Package Size: #59072-1 50ul #59072-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	IRF6 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	IRF6 Antibody detects endogenous levels of IRF6
Immunogen Description	A synthesized peptide derived from human IRF6
Other Names	Interferon regulatory factor 6; IRF 6; LPS;OFC 6; PIT; PPS; PPS1; VWS; VWS1;
Accession No.	Uniprot:O14896
Uniprot	O14896
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

Images



Western blot analysis of IRF6 expression in Jurkat cell lysate.

Product Description

Play an important role in pathogen defense, autoimmunity, lymphocyte development, cell growth, and susceptibility to transformation. Plays a role in regulating mammary epithelial cell proliferation (By similarity). May regulate WDR65 transcription.

Background

Play an important role in pathogen defense, autoimmunity, lymphocyte development, cell growth, and susceptibility to transformation. Plays a role in regulating mammary epithelial cell proliferation (By similarity). May regulate WDR65 transcription.

Note: This product is for in vitro research use only