

Hsp20 Rabbit mAb

Catalog No: #59171

Package Size: #59171-1 50ul #59171-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Hsp20 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC FC
Species Reactivity	Human Mouse Rat
Specificity	Hsp20 Antibody detects endogenous levels of total Hsp20
Immunogen Description	A synthesized peptide derived from human Hsp20
Other Names	HSPB6; HEL55; Hsp20; Heat-shock 20 kDa like-protein p20;
Accession No.	Uniprot:O14558
Uniprot	O14558
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:5000~1:10000 IHC 1:50~1:200 FC 1:50

Images



Western blot analysis of Hsp20 expression in mouse skeletal muscle lysate.

Product Description

HSP20 a small heat shock protein of the HSP20 family. Note: This description may include information from UniProtKB.

Background

HSP20 a small heat shock protein of the HSP20 family. Note: This description may include information from UniProtKB.

Note: This product is for in vitro research use only