

Cip4 Rabbit mAb

Catalog No: #60117



Package Size: #60117-1 50ul #60117-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Cip4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse
Specificity	Cip4 Antibody detects endogenous levels of total Cip4
Immunogen Description	A synthesized peptide derived from human Cip4
Other Names	Cdc42 interacting protein 4; CIP4; DCIP4; hSTP; STOT; STP; Thyroid receptor interacting protein 10; trip10;
Accession No.	Uniprot:Q15642
Uniprot	Q15642
Calculated MW	75kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



Western blot analysis of Cip4 expression in HepG2 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

Product Description

Required for translocation of GLUT4 to the plasma membrane in response to insulin signaling (By similarity).

Background

Required for translocation of GLUT4 to the plasma membrane in response to insulin signaling (By similarity).

Note: This product is for in vitro research use only