

MYPN Antibody FITC Conjugated

Catalog No: #C00718F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MYPN Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC,IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 945-995 1320 derived from human MYPN
Conjugates	FITC
Target Name	MYPN
Other Names	MYOP; RCM4; CMH22; CMD1DD; Myopalladin; 145 kDa sarcomeric protein; MYPN
Accession No.	Swiss-Prot#Q86TC9NCBI Gene ID84665
Uniprot	Q86TC9
GeneID	84665;
Excitation Emission	494nm 518nm
Cell Localization	Cytoplasm, Nucleus
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

Component of the sarcomere that tethers together nebulin (skeletal muscle) and nebulin (cardiac muscle) to alpha-actinin, at the Z lines.

Note: This product is for in vitro research use only