

TBX2 Polyclonal Antibody

Catalog No: #49236



Package Size: #49236-1 50ul #49236-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TBX2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse
Specificity	TBX2 Polyclonal Antibody detects endogenous levels of TBX2 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human TBX2.
Target Name	TBX2
Other Names	TBX2;T-box transcription factor TBX2;T-box protein 2
Uniprot	Q13207
GeneID	6909
Calculated MW	74kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500 - 1:2000. IHC 1:100 - 1:300. ELISA: 1:20000.. IF 1:50-200

Images



Western blot analysis of lysates from HeLa cells, using TBX2 Antibody. The lane on the right is blocked with the synthesized peptide.



Immunohistochemical analysis of paraffin-embedded human brain tumor.

Note: This product is for in vitro research use only