

DCP ANTIBODY FITC Conjugated

Catalog No: #C00726F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DCP ANTIBODY FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt,IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic phosphopeptide aa 60-110 622 derived from human DCP
Conjugates	FITC
Target Name	DCP
Other Names	PT; THPH1; RPRGL2; Prothrombin; Coagulation factor II; F2
Accession No.	Swiss-Prot#P00734NCBI Gene ID2147
Uniprot	P00734
GeneID	2147;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular, Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Background

Thrombin, which cleaves bonds after Arg and Lys, converts fibrinogen to fibrin and activates factors V, VII, VIII, XIII, and, in complex with thrombomodulin, protein C. Functions in blood homeostasis, inflammation and wound healing.

Note: This product is for in vitro research use only