

APC anti-human CD274

Catalog No: #FC0105

Package Size: #FC0105-1 25Tests #FC0105-2 100Tests

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	APC anti-human CD274
Brief Description	mouse monoclonal antibody
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-29E2A3
Isotype	mouse IgG2b, κ
Applications	FC
Species Reactivity	human
Conjugates	APC
Target Name	CD274
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 μ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1×10^6 cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



PHA-stimulated (3 days) human peripheral blood lymphocytes stained with APC anti-human CD274 (colored histogram) or APC mouse IgG2b κ (gray histogram)

Note: This product is for in vitro research use only