

factor V Antibody FITC Conjugated

Catalog No: #C00745F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | factor V Antibody FITC Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Purified by Protein A. |
| Applications | IF |
| Species Reactivity | Hu Ms Rt |
| Immunogen Description | KLH conjugated synthetic peptide aa 300-350 2224 derived from human factor V |
| Conjugates | FITC |
| Target Name | factor V |
| Other Names | FVL; PCCF; THPH2; RPRGL1; Coagulation factor V; Activated protein C cofactor; Proaccelerin, labile factor; F5 |
| Accession No. | Swiss-Prot#P12259NCBI Gene ID2153 |
| Uniprot | P12259 |
| GeneID | 2153; |
| Excitation Emission | 494nm 518nm |
| Cell Localization | Secreted |
| Concentration | 1mg ml |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

IF=1:50-200

Background

Central regulator of hemostasis. It serves as a critical cofactor for the prothrombinase activity of factor Xa that results in the activation of prothrombin to thrombin.

Note: This product is for in vitro research use only