

factor V Antibody HRP Conjugated

Catalog No: #C00745H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	factor V Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB,,IHC-P,IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 300-350 2224 derived from human factor V
Conjugates	HRP
Target Name	factor V
Other Names	FVL; PCCF; THPH2; RPRGL1; Coagulation factor V; Activated protein C cofactor; Proaccelerin, labile factor; F5
Accession No.	Swiss-Prot#P12259NCBI Gene ID2153
Uniprot	P12259
GeneID	2153;
Excitation Emission	N A
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Central regulator of hemostasis. It serves as a critical cofactor for the prothrombinase activity of factor Xa that results in the activation of prothrombin to thrombin.

Note: This product is for in vitro research use only