

FITC anti-mouse H-2Kd/H-2Dd

Catalog No: #FC1035



Package Size: #FC1035-1 25Tests #FC1035-2 100Tests

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description


Product Name	FITC anti-mouse H-2Kd/H-2Dd
Brief Description	mouse monoclonal antibody
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-34-1-2S
Isotype	mouse IgG2a
Applications	FC
Species Reactivity	mouse
Conjugates	FITC
Target Name	H-2Kd/H-2Dd
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Note: This product is for in vitro research use only