



Purified anti-human HLA-DP (MHC Class II)  
Monomorphic

Catalog No: #FC1552

Package Size: #FC1552-1 100ug #FC1552-2 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

Product Name	Purified anti-human HLA-DP (MHC Class II) Monomorphic
Brief Description	mouse monoclonal antibody
Host Species	mouse
Clone No.	s-B7/21
Species Reactivity	human
Storage	Store at 4°C.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only