

Purified anti-mouse CD278

Catalog No: #FC1643



Package Size: #FC1643-1 100ug #FC1643-2 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---------------------------|
| Product Name | Purified anti-mouse CD278 |
| Brief Description | rat monoclonal antibody |
| Host Species | rat |
| Clone No. | s-7E.17G9 |
| Species Reactivity | mouse |
| Storage | Store at 4°C. |

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only