

Purified anti-mouse MAdCAM-1

Catalog No: #FC1671



Package Size: #FC1671-1 100ug #FC1671-2 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Purified anti-mouse MAdCAM-1
Brief Description	rat monoclonal antibody
Host Species	rat
Clone No.	s-MECA-367
Species Reactivity	mouse
Storage	Store at 4°C.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only