

## NSUN3 Rabbit pAb

Catalog No: #62025



Package Size: #62025-1 50ul #62025-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	NSUN3 Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	affinity purified by Protein A
Applications	IHC-P IHC-F ICC/IF
Species Reactivity	Human, Mouse
Immunogen Description	KLH conjugated synthetic peptide derived from human NSUN3
Other Names	6720484A09Rik; AU022521; FLJ22109; FLJ22609; MST077; MSTP077; NOL1/NOP2/Sun domain family 3; NOL1/NOP2/Sun domain family member 3; NOP2/Sun domain family member 3; NSUN3; NSUN3_HUMAN; OTTMUSP00000028633; Putative methyltransferase NSUN3; UG0651E06.
Uniprot	Q9H649
GeneID	63899
Calculated MW	38 kDa
Concentration	1 mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-P=1:100-500, IHC-F=1:100-500, ICC/IF=1:100-500

Note: This product is for in vitro research use only