

Anti-PINK1 (Phospho-S228) Antibody

Catalog No: #62030



Package Size: #62030-1 50ul #62030-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Anti-PINK1 (Phospho-S228) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB, IH
Species Reactivity	H, M, R
Immunogen Description	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S228 of human PINK1 protein.
Other Names	Serine/threonine-protein kinase PINK1 mitochondrial; BRPK; PTEN-induced putative kinase protein 1
Uniprot	Q9BXM7
GeneID	65018
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	store at -20°C for one year

Application Details

WB (1/500 - 1/1000), IH (1/50 - 1/100)

Note: This product is for in vitro research use only