

## COX14 rabbit pAb

Catalog No: #62033



Package Size: #62033-1 50ul #62033-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	COX14 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of human COX14
Immunogen Description	Synthesized peptide derived from human COX14
Target Name	COX14 C12orf62
Uniprot	Q96I36
GeneID	84987
Calculated MW	6kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only