

COQ2 rabbit pAb

Catalog No: #62061



Package Size: #62061-1 50ul #62061-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	COQ2 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	This antibody detects endogenous levels of COQ2 at Human/Mouse/Rat
Immunogen Description	Synthesized peptide derived from human COQ2
Target Name	COQ2 CL640
Uniprot	Q96H96
GeneID	27235
Calculated MW	41kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000

Images



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4°over night

Note: This product is for in vitro research use only