

## DCAF8 Rabbit pAb

Catalog No: #62065



Package Size: #62065-1 50ul #62065-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	DCAF8 Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	Human Mouse Rat
Species Reactivity	WB
Specificity	This antibody detects endogenous levels of human DCAF8
Immunogen Description	Synthesized peptide derived from human DCAF8
Target Name	DCAF8 H326 WDR42A
Uniprot	Q5TAQ9
GeneID	50717
Calculated MW	66kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only