

## DAPL1 Rabbit pAb

Catalog No: #62069

Package Size: #62069-1 50ul #62069-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	DAPL1 Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC;IF
Species Reactivity	Human;Mouse
Specificity	This antibody detects endogenous levels of DAPL1 at Human/Mouse
Immunogen Description	Synthesized peptide derived from human DAPL1
Target Name	DAPL1 EEDA
Uniprot	A0PJW8
GeneID	92196
Calculated MW	12kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

IHC 1:50-200; IF 1:50-200

## Images



Immunohistochemical analysis of paraffin-embedded human uterus.

Note: This product is for in vitro research use only