

CTHRC1 Rabbit pAb

Catalog No: #62082



Package Size: #62082-1 50ul #62082-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | CTHRC1 Rabbit pAb |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications | WB |
| Species Reactivity | Human;Mouse;Rat |
| Specificity | This antibody detects endogenous levels of CTHRC1 at Human, Mouse,Rat |
| Immunogen Description | Synthesized peptide derived from human CTHRC1 |
| Target Name | CTHRC1 UNQ762/PRO1550 |
| Uniprot | Q96CG8 |
| GeneID | 115908 |
| Calculated MW | 27kD |
| Concentration | 1 mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C/1 year |

Application Details

WB 1:500-2000

Note: This product is for in vitro research use only