

Recombinant Human CPEB3/CPE-BP3 Protein, N-His

Catalog No: #AP98555



Package Size: #AP98555-1 50ug #AP98555-2 100ug #AP98555-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CPEB3/CPE-BP3 Protein, N-His
Brief Description	Recombinant Protein
Host Species	E.coli
Other Names	Cytoplasmic polyadenylation element-binding protein 3, CPE-BP3, CPE-binding protein 3, hCPEB-3, CPEB3
Uniprot	Q8NE35
Target Species	Homo sapiens (Human)
SDS-PAGE MW	28.61 kDa
Tag Info	N-His
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

>90% as determined by SDS-PAGE.

Note: This product is for in vitro research use only