

Collagen VIII alpha 1 Rabbit pAb

Catalog No: #62072

Package Size: #62072-1 50ul #62072-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Collagen VIII alpha 1 Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB;IHC-P;IHC-F;ICC/IF
Species Reactivity	Human;Mouse;Rat
Immunogen Description	KLH conjugated synthetic peptide derived from human Collagen VIII alpha 1
Other Names	C3orf7; CO8A1_HUMAN; COL8A1; Collagen alpha 1(VIII) chain [Precursor]; Endothelial collagen; MGC9568; Vastatin.
GeneID	1295
Calculated MW	19/79kDa
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500

Note: This product is for in vitro research use only