

Monkey Apolipoprotein E (APOE) ELISA Kit

Catalog No: #EK20868

Package Size: #EK20868-1 48T #EK20868-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Monkey Apolipoprotein E (APOE) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Monkey
Other Names	AD2;LDLCQ5; LPG;MGC1571;apolipoprotein E3

Application Details

Detect Range:7.8 ng/mL - 500 ng/mL

Sensitivity:3.9 ng/mL

Sample Type:Serum, Plasma, Other biological fluids

Sample Volume: 1-200 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only