

GLUT-1 Mouse mAb

Catalog No: #62074

Package Size: #62074-1 50ul #62074-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GLUT-1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a,Kappa
Purification	The antibody was affinity-purified from ascites by affinity-chromatography using specific immunogen.
Applications	IHC;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	The antibody can specifically recognize human GLUT-1 protein.
Immunogen Description	Synthesized peptide derived from human GLUT-1
Target Name	SLC2A1
Other Names	GLUT1
Uniprot	P11166
GeneID	6513
Calculated MW	54kD
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage	Store at -20°C/1 year

Application Details

IHC 1:200-400; IF 1:50-200; ELISA 1:500-5000

Images



Mouse kidney tissue was stained with Anti-GLUT-1 Antibody

Note: This product is for in vitro research use only