

Caspase-1 Mouse mAb

Catalog No: #62079



Package Size: #62079-1 50ul #62079-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Caspase-1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b,Kappa
Purification	Protein G
Applications	WB;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of Caspase-1 protein.
Immunogen Description	Synthesized peptide derived from human Caspase-1
Target Name	CASP1 IL1BC IL1BCE
Uniprot	P29466
GeneID	834
Calculated MW	45kDa
Concentration	1 mg/ml
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000; IF 1:100-500; ELISA 1:1000-5000

Images



Whole cell lysates were separated by 10% SDS-PAGE, and the membrane was blotted with anti-Caspase-1 antibody.

Note: This product is for in vitro research use only