

GSDMD Mouse mAb

Catalog No: #62096



Package Size: #62096-1 50ul #62096-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


Description

Product Name	GSDMD Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1,kappa
Purification	The antibody was affinity-purified from ascites by affinity-chromatography using specific immunogen.
Applications	WB;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of GSDMD protein.
Target Name	GSDMD
Uniprot	P57764
GenelD	53kD,31kD
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000; IF 1:100-500; ELISA 1:1000-5000

Images



Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-GSDMD antibody.

Note: This product is for in vitro research use only