

Recombinant Human ZNF207 Protein, N-GST & C-His

Catalog No: #AP98556



Package Size: #AP98556-1 50ug #AP98556-2 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ZNF207 Protein, N-GST & C-His
Brief Description	Recombinant Protein
Other Names	BUB3-interacting and GLEBS motif-containing protein ZNF207, BuGZ, hBuGZ, Zinc finger protein 207, ZNF207
Uniprot	O43670
Target Species	Homo sapiens (Human)
SDS-PAGE MW	39.02 kDa
Tag Info	N-GST & C-His
Target Sequence	Asp21-Ser113
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

>90% as determined by SDS-PAGE.

Note: This product is for in vitro research use only