

# Mouse Transcription factor SOX-9 (SOX9) ELISA Kit

Catalog No: #EK20892



Package Size: #EK20892-1 48T #EK20892-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                   |
|--------------------|---------------------------------------------------|
| Product Name       | Mouse Transcription factor SOX-9 (SOX9) ELISA Kit |
| Brief Description  | ELISA Kit                                         |
| Applications       | ELISA                                             |
| Species Reactivity | Mouse                                             |

## Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.103ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only