

IL-17 Rabbit pAb

Catalog No: #62241

Package Size: #62241-1 50ul #62241-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-17 Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;Flow Cyt;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of Human IL-17
Immunogen Description	Synthesized peptide derived from human IL-17
Target Name	IL17A CTLA8 IL17
Other Names	Interleukin-17A ; IL-17 ; IL-17A ; Cytotoxic T-lymphocyte-associated antigen 8 ; CTLA-8
Uniprot	Q16552
GeneID	3605
Calculated MW	17kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000; Flow Cyt 1:50-200; ELISA 1:5000-20000

Note: This product is for in vitro research use only