

Recombinant Human FGF-18 (C-6His)

Catalog No: #AP89690




Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Recombinant Human FGF-18 (C-6His)
Brief Description	Recombinant Protein
Host Species	Mammalian Cell
Uniprot	O76093
Calculated MW	25-37kDa, reducing conditions
Tag Info	C-6His
Target Sequence	Glu28-Ala207
Formulation	Lyophilized from a 0.2 M filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at β -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at β -20°C for 3 months.

Images



Greater than 95% as determined by reducing SDS-PAGE.



Immobilized Recombinant Human FGF-18 at 5ug/ml (100 ul/well) can bind Recombinant Rat FGFR4. The ED50 of Recombinant Rat FGFR4 is 0.196 ug/ml. (Regularly tested).

Note: This product is for in vitro research use only