

DHX37 Antibody

Catalog No: #62175

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DHX37 Antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total DHX37 protein
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human DHX37
Target Name	DHX37; DDX37; KIAA1517;
Uniprot	Q8IY37;
Calculated MW	Predicted band size: 127 kDa
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-2000

Note: This product is for in vitro research use only