

## SOSB1 Antibody

Catalog No: #62177

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                                                       |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name          | SOSB1 Antibody                                                                                                        |
| Brief Description     | Rabbit Polyclonal                                                                                                     |
| Host Species          | Rabbit                                                                                                                |
| Clonality             | Polyclonal                                                                                                            |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications          | WB                                                                                                                    |
| Species Reactivity    | Human; Mouse; Rat                                                                                                     |
| Specificity           | The antibody detects endogenous level of total SOSB1 protein                                                          |
| Immunogen Type        | Peptide                                                                                                               |
| Immunogen Description | Synthesized peptide derived from human SOSB1                                                                          |
| Target Name           | NABP2; OBFC2B; SSB1; LP3587;                                                                                          |
| Uniprot               | Q9BQ15;                                                                                                               |
| Calculated MW         | Predicted band size: 23 kDa                                                                                           |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| Storage               | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.                                       |

## Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only